
GASEX OEM HD-05-xx

Arcoptix GASEX OEM HD-05 产品简介:

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

FTIR	
Design	Permanently aligned, maintenance free, hermetic.
Resolution [cm⁻¹]	0.5
Beamsplitter / window material	ZnSe
Spectral range [cm⁻¹]	5000 - 830 (5000 – 1800 - optional)
Detector	MCT (4-TE cooled)
Detectivity D* [cm Hz^{1/2} W⁻¹]	> 2.5x10 ⁹ (4x10 ¹⁰ - optional)
Light source	Broadband SiC, 1550 K
Reference laser	Stabilized semiconductor laser
Scan frequency [s⁻¹]	0.7
Signal-to noise ratio	> 3,500:1 (for one single scan)
Mounting orientation (operation)	Position independent
Absotrans^(TM)	Active suppression of H ₂ O and CO ₂
GAS CELL	
Path length	5m
Internal volume	0.2L
Transmission	>50%
Temperature range [°C]	20 to 200 (not condensing)
Construction	Aluminum with inert coating.
Mirrors parabolic,	Rhodium protected, gold coated aluminum
Windows	BBAR ZnSe
Gas connectors	For 6mm (Swagelok) or ¼” diameter tubing (customizable)
SYSTEM	
Communication	USB 3.0 or Ethernet
Operating system	Windows 7/8/10 and Linux
API	Software or hardware server for data retrieval, GUI for spectral acquisition, Optionally Panorama by Labcognition.
Power consumption	FTIR only: 32W @ 12VDC (adapter supplied) Gas cell heaters: 400W (peak) 20W (in the heated state) @ 110-230VAC (50-60Hz) or 12 V (customizable)
Dimensions [mm]	380 x 180 x 120
Weight	3.9 kg

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Examples of CO spectrum measured with GASEX HD-05-12 in comparison to the GASEX ST-20-12 with 2 cm⁻¹ (inset of the picture) resolution.


案例2:用GASEX HD-05测量的120ppm 氯化氢 (Hydrogen Chloride) 光谱

Examples of HCl measurements using a GASEX HD-05-12 together with HITRAN simulation of the spectra for HCl with ^{35}Cl and ^{37}Cl isotopes in the molecule. The isotope ratio can be determined.

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FTIR – gas cell unit for gas analyzers

FTIR – gas cell unit

Examples of various gases spectra measured using using a GASEX HD05-12 (screenshot) showing a variety of substances measured, within the 2-12 microns range or operation of these spectrometers.

Substance	LOD* [ppm]	Substance	LOD* [ppm]
HF	0.18	HBr	0.18
HCl	0.14	CO2	0.15
HCN	0.32	n-Octane	0.14
NH3	0.12	Formaldehyde	0.21
CO	0.2	Acetone	0.28
N2O	0.03	Phosgene	0.013
NO	1**	Chloroform	0.032
NO2	0.26	Methanol	0.43
SO2	0.03	Ethanol	0.14

* Values estimated from as three times standard deviation for 64 spectra averaged

** Measured with one single scan using a short wavelength - 6 micron cut-off detector.

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